OPERATION MANUAL

LLG-uniSTIRRER 7 pro

Thank you very much for purchasing this LLG-uniSTIRRER 7 pro. Please read the following instructions carefully and always keep this manual within easy reach.

Contents

1. Preface/Service	3	7. Working with external temperature sensor	11
2. Safety Instructions	3	8. Residual heat warning (HOT)	12
3. Proper use	5	9. Faults	12
4. Inspection	6	10. Maintenance and Cleaning	13
4.1 Control elements	6	11. Associated standards and regulations	14
4.2 Display	6	12. Specifications	15
5. Control	7	13. Ordering information	16
5.1 Steuerelemente	7	14. Warranty	17
5.2 Display	9		
6. Trial run	10		

1. Preface/Service

Thank you very much for purchasing the LLG-uniSTIRRER 7 pro. You should read this Manual carefully, follow the instructions and procedures, and be aware of all the cautions when using this instrument.

Please provide the customer care representative with the following information:

- Serial number (on the rear panel)
- Certification
- Description of problem (i.e., hardware or software)
- Methods and procedures adopted to resolve the problems
- Your contact information

2. Safety Instructions

	<p>Warning!</p> <ul style="list-style-type: none"> • Read the operating instructions carefully before use. • Ensure that only trained staff works with the instrument.
	<p>Risk of burn!</p> <ul style="list-style-type: none"> • Caution when touch the housing parts and the hotplate (heating model) which can reach temperature of 550 °C. • Pay attention to the residual heat after switching off.
	<p>Protective ground contact!</p> <ul style="list-style-type: none"> • Make sure that socket must be grounded (protective ground contact) before use.

- When working wear personal safety guards to avoid the risk from:
 - Splashing and evaporation of liquids
 - Release of toxic or combustible gases
- Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- Gradually increase the speed, reduce the speed if:
 - Stirring bar breaks away due to high speed
 - The instrument is not running smoothly, or container moves on the base plate
- Temperature must always be set to at least 50 °C lower than the fire point of the media used.
- Be aware of hazards due to:
 - Flammable materials or media with a low boiling temperature
 - Overfilling of media
 - Unsafe container
- Process pathogenic materials only in closed vessels.
- Check the instrument and accessories prior to each use. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter. Accessories must be securely attached to the device and can not come off by themselves. Always disconnect the plug before fitting accessories.
- When the external temperature sensor needed, the tip of the measuring sensor must be at least 5-10 mm from vessel bottom and wall.
- The instrument can only be disconnected from the main power supply by pulling out the main or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the main power supply cable does not touch the hotplate. Do not cover the device.
- Forbid to put pressure and over heat media on the surface of glass ceramic that can be caused surface broken.
- The instrument may only be opened by experts.
- Keep away from high magnetic field.

3. Proper use

The LLG-uniSTIRRER 7 pro is designed for mixing and/or heating liquids in schools, laboratories or factories.

- Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 100 mm)


Figure 1

This device is not suitable for using in residential areas or other constraints mentioned in Chapter 1.

4. Inspection

4.1 Receiving Inspection

Unpack the equipment carefully and check for any damages which may have arisen during transport. Please contact your supplier for technical support.

4.2 Listing of Items

The package includes the following items:

Items	Qty
Main unit	1
Power cable	1
User Manual	1
Temperature sensor	1
Support rod	1
Clamp for Temperature sensor	1

Note:

If there is any apparent damage to the system, please do not plug it into the power line.

5. Control

5.1 Control elements


Figure 2
LLG-uniSTIRRER 7 pro


Figure 3
Power Switch

Items	Descriptions
Speed control knob Stir	Set the rotary speed. The stirring function is switched ON or OFF by pushing the knob.
Temperature control knob Heat (Hotplate)	Set the temperature. The heating function is switched ON or OFF by pushing the knob.
LCD Display	LCD displays the real working state and all settings.
LED Heat (Hotplate)	When the heating function is switched ON, the LED Heat is lit.
LED Stir	When the stirring function is switched ON, the LED Stir is lit.
Power Switch	Switch ON or OFF the instrument.

5.2 Display


Figure 4
LCD digital model

Characters	Descriptions
Temp und °C	Display temperature when the heating function is switched ON.
Stir	Display stirring state when the stirring function is switched ON.
Hot	Display hot warning if the heating plate temperature is above 50 °C after switching OFF the heating function.
Probe	Display when using external probe.
Remote	Display in case of remote control.
Err	Display in case of remote control.
Rated value / Real value	Display value in case of heating or stirring function switching ON.

Note!

If both heating and stirring functions have been started at the same time, heating function always has higher priority. If in this case speed is changed via the stirring knob, it displays stirring speed and reverses to temperature in the duration of 5 second

6. Trial run

- Make sure the required operating voltage and power supply voltage match.
- Ensure the socket must be properly grounded.
- Plug in the power cable, ensure the power is on and begin initializing.
- Add the medium into the vessel with an appropriate stirring bar.
- Place vessel on the work plate.
- Set the target stirring speed and begin.
- Set the target temperature and start heating (hotplate model).
- Stop the heating and stirring functions.

If these operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact your supplier for technical support.

Warning!

Forbid to transfer the vessel when the instrument working.

7. Working with external temperature sensor

The external temperature sensor PT1000 is manufacture's standard accessory. If the sensor is plugged in, "Probe" will be shown on the LCD digital display to indicate the sensor is operating. The setting value of external temperature sensor and actual temperature are displayed. Safe circuit controls hotplate temperature. Comparing with the temperature control of the hotplate, the external temperature sensor can control the medium's temperature more precise. The heating function will be stopped automatically under abnormal conditions. Please operate follow the instructions below:

- Switch OFF the instrument.
 - Ensure the external temperature sensor is inserted in the media heated.
 - Switch ON the instrument and run heating function.
- If the heating function does not work, please contact your supplier for technical support.


Figure 5

8. Residual heat warning (HOT)

In order to prevent the risk of burns from the hotplate, digital hotplate model has a residual heat warning function. When the heating function is switched off and the heating plate temperature is still above 50 °C, "Hot" will flash to warn that there is a hazard of burns from the hotplate. When the hotplate temperature drops to below 50 °C, the unit will automatically switch off. If users want to turn off the LCD or LED screen immediately, just pull out the plug directly. When the power OFF, the residual heat warning function cannot be run.

9. Faults

- Instruments can't be powered ON
 - Check whether the power line is unplugged
 - Check whether the fuse is broken or loose
- Fault in power ON self test
 - Switch OFF the unit, then switch ON and reset the instruments to factory default setting.
- Stir speed cannot reach set point
 - Excessive medium viscosity may cause abnormal speed reduction of the motor
- Unit cannot be powered OFF when switched OFF.
 - Check if the residual heat warning function is still ON and hotplate temperature is above 50 °C (the LCD/LED screen still work and "Hot" flash).

If these faults are not resolved, please contact the supplier

10. Maintenance and Cleaning

- Proper maintenance can keep instruments working properly and lengthen its lifetime.
- Do not spray cleanser into the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use recommended cleansers:

Dyes	Isopropyl alcohol
Construction materials	Water containing tenside / Isopropyl alcohol
Cosmetics	Water containing tenside / Isopropyl alcohol
Foodstuffs	Water containing tenside
Fuels	Water containing tenside

Wear the proper protective gloves during cleaning of the instrument.

- Before using other method for cleaning or decontamination, the user must ascertain with your supplier that this method will not harm the instrument.
- Send in the case of service the instrument back in the packaging carton. Storage packing is not sufficient for the back dispatch. Use additionally a suitable transportation packing.
- The enamel makes the hotplate easier to care for and more resistant to acids and bases. Because of it, however, the hotplate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

Warning!

Power OFF when maintenance and cleaning.

11. Associated standards and regulations

Construction in accordance with the following safety standards:	EN 61010-1 UL 3101-1 CAN/CSA C22.2(1010-1) EN 61010-2-10
Construction in accordance with the following EMC standards:	EN 61326-1
Associated EU guidelines:	EMC-guidelines: 89/336/EWG Instrument guidelines: 73/023/EWG

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

12. Specifications

Items	Specifications LLG-uniSTIRRER 7 pro
Voltage [VAC]	100-120/200-240
Frequency [Hz]	50/60
Power [W]	1050
Stirring point position quantity	1
Max. stirring quantity (H ₂ O) [l]	20
Max. magnetic bar [L×Ø, mm]	80×10
Motor type	DC brushless motor
Max. power input of motor [W]	18
Max. power output of motor [W]	10
Speed range [rpm]	100-1500
Rotary speed display	LCD
Plate material	Glass ceramic
Dimensions of workplate (mm)	184×184
*Heating power [W]	1000
*Temperature range [°C]	RT-550, increment: 1

Items	Specifications LLG-uniSTIRRER 7 pro
Temperature display accuracy [°C]	±0.1
Control accuracy of heating temperature [°C]	±1 (100 °C below) ±1 % (100 °C above)
The safety temperature of the hotplate [°C]	580
Temperature sensor in medium	PT1000
Control accuracy of heating temperature with temperature sensor [°C]	±0.2
Residual heat warning	50 °C
Dimensions (mm)	215×360×112
Weight [kg]	5.3
Permitted ambient temperature [°C]	5-40
Permitted relative humidity	80 %
Protection class acc. to DIN 60529	IP21
RS232 interface	Yes

13. Ordering information

Cat No.	Descriptions
6.263 448	LLG-uni <i>STIRRER</i> 7 pro with EU plug
6.263 449	LLG-uni <i>STIRRER</i> 7 pro with UK plug
6.263 415	Temperature sensor for LLG-uni <i>STIRRER</i> 7 pro, length 230 mm, sensor length 198 mm
6.263 417	Temperature sensor glass coated for LLG-uni <i>STIRRER</i> 7 pro, length 230 mm, sensor length 198 mm
6.263 416	Rod and support clamp of PT1000

14. Warranty

This product is warranted to be free from defects in material and workmanship for a period of three (3) years from date of purchase.

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in this instruction manual. This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship. This warranty does not cover any incidental or consequential damages, commercial loss or any other damages from the use of this product.

The warranty does not cover damage to paint or finish and defects or damages caused by physical and chemical abuse or normal wear and tear. The warranty is invalidated by any nonfactory modification, which will immediately terminate all liabilities on us for the products or damages caused by its use. The buyer and its customer shall be responsible for the product or use of products as well as any supervision required for safety. If requested the products must be returned to the distributor in well packed and insured manner and all shipping charges must be paid.

Information on the disposal of electrical and electronic devices in the European Community:

The disposal of electrical devices is regulated within the European Community by national regulations based on the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). According to these regulations, any devices supplied after 13.06.2005 in the business to business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste. They are marked with the following symbol to indicate this.

As disposal regulations within the EU may vary from country to country, please contact your supplier if necessary.

PRODUCT DISPOSAL

In case the product is to be disposed of, the relevant legal regulations are to be followed.

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Lab Logistics Group GmbH