

**Liner di ingresso per GC Perkin-Elmer****Trajan Scientific**

Liner di ingresso per GC Perkin-Elmer. Liner per Perkin-Elmer Autosystem e Clarus 500, 600. I liner di ingresso SGE vengono forniti come una unità completa, confezionata singolarmente. 5 o 25 confezioni, individualmente confezionate. Da completare con o-ring ed anelli di tenuta appropriati. Ogni confezione viene fornita con il risultato dei test di qualità.



| Code        | Description                                                                | Packaging |
|-------------|----------------------------------------------------------------------------|-----------|
| LLG06206049 | Split/Splitless Focus Liner, Øest.6,2mm, Øint.4mm, lungh. 92 mm            | 5 pz      |
| LLG06206050 | Split/Splitless Tapered Focus Liner, Øest.6,2mm, Øint.4mm, lungh.92mm      | 5 pz      |
| LLG07670774 | Liner AG 4mm ID, GN QW                                                     | 25 pz     |
| LLG09003567 | Split straight through liner, Øest.6,2mm, Øint.4mm, lungh.92mm             | 5 pz      |
| LLG09003568 | Split/Splitless Focus Liner per PSS Injector, Øest.4mm, Øint.2mm, L.86,2mm | 5 pz      |
| LLG09003569 | PTV liner con restrizione ID 0,25mm, Øest.2mm, Øint.1mm, lungh.88mm        | 5 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

