




Colonne GC per ammine

Macherey Nagel

Colonne GC per ammine



Code	Description	Packaging
MNL73311025	FS capillary column CW 20 M-AM ID: 0.25 mm, film thickness: 0.25 µm length: 25 m	1 pz
MNL73311050	FS capillary column CW 20 M-AM ID: 0.25 mm, film thickness: 0.25 µm length: 50 m	1 pz
MNL73311110	FS capillary column CW 20 M-AM ID: 0.10 mm, film thickness: 0.25 µm length: 10 m	1 pz
MNL73329925	FS capillary column CW 20 M-AM ID: 0.32 mm, film thickness: 0.25 µm length: 25 m	1 pz
MNL73329950	FS capillary column CW 20 M-AM ID: 0.32 mm, film thickness: 0.25 µm length: 50 m	1 pz
MNL73344250	FS-Capillary column CW 20 M-AM ID: 0.32 mm, film thickness: 0.35 µm length: 50 m	1 pz
MNL73355125	FS capillary column CW 20 M-AM ID: 0.53 mm, film thickness: 1.0 µm length: 25 m	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

