



Optima FFAP plus colonne capillari polari per GC

Macherey Nagel

Optima FFAP plus colonne capillari polari per GC



Code	Description	Packaging
MNL72624130	Capillary column OPTIMA FFAPplus ID: 0.25 mm, film thickness: 0.25 µm length: 30 m pack of 1	1 pz
MNL72624160	Capillary column OPTIMA FFAPplus ID: 0.25 mm, film thickness: 0.25 µm length: 60 m pack of 1	1 pz
MNL72624230	Capillary column OPTIMA FFAPplus ID: 0.25 mm, film thickness: 0.50 µm length: 30 m pack of 1	1 pz
MNL72624260	Capillary column OPTIMA FFAPplus ID: 0.25 mm, film thickness: 0.50 µm length: 60 m pack of 1	1 pz
MNL72624330	Capillary column OPTIMA FFAPplus ID: 0.32 mm, film thickness: 0.25 µm length: 30 m pack of 1	1 pz
MNL72624360	Capillary column OPTIMA FFAPplus ID: 0.32 mm, film thickness: 0.25 µm length: 60 m pack of 1	1 pz
MNL72624630	Capillary column OPTIMA FFAPplus ID: 0.32 mm, film thickness: 0.50 µm length: 30 m pack of 1	1 pz
MNL72624660	Capillary column OPTIMA FFAPplus ID: 0.32 mm, film thickness: 0.50 µm length: 60 m pack of 1	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

