



SPE glass columns, CHROMABOND Florisil, 150–250 μm , 6 mL/2000 mg

MACHERY NAGEL

Florisil is used for sample preparation (solid phase extraction) of e.g., organic tin compounds and aliphatic carboxylic acids. Florisil is a magnesium silicate.



Code	Description	Packaging
MNL730239G	Chromab. columns Florisil, 6 mL, 2000 mg	30 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

