



## Acido nitrico 65%

Synonyms:

### Classification

### Danger

Molecular weight **63.013** ONU**2031**CAS **7697-37-2** Transport hazard class:**8**

EEC-N: Packing

**II****Acido nitrico 65%****> RPE - Per analisi - ISO****Grado Analitico**

|                      |                         |     |
|----------------------|-------------------------|-----|
| Description          | Clear colourless liquid | -   |
| Identification       | Positive                | -   |
| Density d20/4        | 1.380 - 1.400           | -   |
| Chloride             | $\leq 0.1$              | ppm |
| Phosphate            | $\leq 0.5$              | ppm |
| Heavy metals (as Pb) | $\leq 0.2$              | ppm |



|                     |          |     |
|---------------------|----------|-----|
| Residue on ignition | <= 4     | ppm |
| Silicate            | <= 1     | ppm |
| Sulfate             | <= 0.5   | ppm |
| Ag                  | <= 0.02  | ppm |
| Al                  | <= 0.05  | ppm |
| As                  | <= 0.005 | ppm |
| Ba                  | <= 0.1   | ppm |
| Be                  | <= 0.02  | ppm |
| Bi                  | <= 0.1   | ppm |
| Ca                  | <= 0.5   | ppm |
| Cd                  | <= 0.005 | ppm |
| Co                  | <= 0.01  | ppm |
| Cr                  | <= 0.1   | ppm |
| Cu                  | <= 0.01  | ppm |
| Fe                  | <= 0.2   | ppm |
| K                   | <= 0.1   | ppm |
| Li                  | <= 0.02  | ppm |
| Mg                  | <= 0.1   | ppm |
| Mn                  | <= 0.01  | ppm |
| Mo                  | <= 0.05  | ppm |



|                     |         |     |
|---------------------|---------|-----|
| Na                  | <= 0.5  | ppm |
| Ni                  | <= 0.05 | ppm |
| Pb                  | <= 0.02 | ppm |
| Sr                  | <= 0.02 | ppm |
| Ti                  | <= 0.1  | ppm |
| Tl                  | <= 0.05 | ppm |
| V                   | <= 0.05 | ppm |
| Zn                  | <= 0.05 | ppm |
| Assay (acidimetric) | 64 - 66 | %   |

| ID                 | Size      | Packaging                           |
|--------------------|-----------|-------------------------------------|
| 000000000000408021 | 1 l       | Bottiglia di vetro rivestita di PVC |
| 408021X6           | 6 x 1 L   | Bottiglia di vetro rivestita di PVC |
| 000000000000408022 | 1 l       | Bottiglia di vetro                  |
| 408022X6           | 6 x 1 L   | Bottiglia di vetro                  |
| 000000000000408025 | 2.5 l     | Bottiglia di vetro                  |
| 408025X4           | 4 x 2,5 L | Bottiglia di vetro                  |
| 000000000000408027 | 34 kg     | Tanica di plastica                  |
| 000000000000524535 | 1 l       | Bottiglia di plastica               |
| 524535X6           | 6 x 1 L   | Bottiglia di plastica               |



| Acido nitrico 65% > RPE - Per analisi - ISO |           | Grado Analitico       |
|---------------------------------------------|-----------|-----------------------|
| ID                                          | Size      | Packaging             |
| 000000000000524536                          | 2.5 l     | Bottiglia di plastica |
| 524536X4                                    | 4 x 2,5 L | Bottiglia di plastica |



CARLO ERBA Reagents S.r.l.  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA Reagents  
operates with a Certified  
Quality Management  
System

