



Micro capillaries

Micro capillaries

These disposable micro capillaries fill automatically from end to end. Their accuracy is better than 1%. A capillary tube holder is supplied with each pack. Also contains a small bulb with aperture which can be inserted into capillary and used as an effective support for filling and emptying tubes.



Code	Description	Packaging
LLG09020192	Microcapillaries, cap. 0,5 µl, pack o/100	100 pz.
LLG09020193	Microcapillaries, cap. 1,0 µl, pack o/100	100 pz.
LLG09020194	Microcapillaries, cap. 2,0 µl, pack o/100	100 pz.
LLG09020195	Microcapillaries, cap. 5,0 µl, pack o/100	100 pz.
LLG09020196	Microcapillaries, cap. 10,0 µl, pack /100	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

