



Accessori per camera per gel elettroforesi MultiSUB Choice

Thistle Scientific

Accessori per camera per gel elettroforesi MultiSUB Choice.



Code	Description	Packaging
LLG06229543	Vassoio-UV, 150x100 mm	1 pz
LLG06231475	Vassoio UV, 150x150 mm	1 pz
LLG06233618	Pettine tipo MS15-20-1.5, 20 campioni, 1,5 mm, 32 µl	1 pz
LLG06234178	Vassoio UV, 150 x 70 mm	1 pz
LLG06234451	Pettine tipo MS15-20-1, 20 campioni, 1mm, 21 µl	1 pz
LLG06238277	Pettine tipo MS15-35-1, 35 campioni, 1 mm, 10 µl	1 pz
LLG06258162	Pettine per 16 campioni, spessore 1,5 mm, volume campione 41 µL	1 pz
LLG07629621	Bordi per vassoio UV	1 pz
LLG07673222	Pettine 16 campioni MC, Volume Campione 27 µl, Spessore 1,00 mm	1 pz
LLG09584666	Pettine tipo MS15-30MC-1.5, 30 campioni MC, 1,5 mm, 19 µl	1 pz
LLG09584815	Pettine 28 campioni MCSS, volume campione 11 µl, spessore 1,00 mm, per vasca elettroforesi MultiSUB Choice	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

