

**SPE columns, CHROMABOND HLB, 30 μ m, 3 mL/200 mg****MACHERY NAGEL**

State-of-the-art spherical polymer for highest demands in solid phase extraction. Hydrophilic-lipophilic balanced (HLB) N-vinylpyrrolidone-divinylbenzene copolymer. Applicable for a wide range of analyte polarities. Enhanced retention for polar compounds. High loadability and outstanding performance. Water wettable – even if bed runs dry, SPE can be continued.



Code	Description	Packaging
MNL730924P30	Chromab. columns HLB (30 µm), 3 mL,200mg	30 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

