

**EC HPLC column (analytical), NUCLEODUR C18 HTec, 3 μ m, 150x2 mm****MACHEREY NAGEL**

Due to high loadability and further characteristics like excellent stability these HPLC columns packed with high purity NUCLEODUR C18 HTec silica gel are the choice for up-scaling to preparative scale.



Code	Description	Packaging
MNL76032520	EC 150/2 NUCLEODUR C18 HTec, 3 µm	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

