

**EC UHPLC column (analytical), NUCLEODUR PolarTec, 1.8 μ m, 30x3 mm****MACHERY NAGEL**

Based on high-purity NUCLEODUR silica gel, a polar group was embedded in a hydrophobic C18 modification. As a result, this HPLC phase shows steric and an increased polar selectivity. This HPLC phase shows a good performance under high aqueous conditions. It is a hydrophobic phase with pronounced hydrophilic properties. Separation mechanism based on hydrophobic (van der Waals) and polar interactions. Suitable for LC/MS due to low bleeding characteristics.



Code	Description	Packaging
MNL76046130	EC 30/3 NUCLEODUR PolarTec, 1.8 µm	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

